







Frequency	Depth	Mode	Applications
3.5-5MHz	90-305mm	B Mode B+M PW Doppler Color Doppler Power Doppler	Abdominal Emergency Cardiac Urology Obstetrics / Gynecology, Lung Kidney

Download our ultrasound software from App Store or Google Play



- √ Unlimited users
- √ Free software
- √ No hidden fees
- √ Support printing reports
- ✓ Chinese, English, Russian, ItalianSpanish, Portuguese (Brazil)

Software Controls

GN(gain)

D(depth)

ENH(enhancement)

DR(dynamic range)

F(frequency)

FocusPos

PRF

WF

Mode

8TGC

Biopsy

Annote

Page left/right up/down

Measurement

B: Length, Area/Circumference,

Angle, Trace, Distance GA (CRL, BPD, GS, FL, HC, AC)

EFW(BPD,FL)

B+M: Heart Rate, Time, Distance

B+PW:Velocity, Heart Rate(2), S/D, Depth



Package includes:

1* wireless probe

1* wireless charger

1* USB charger

1* wrist strap

1* plug

1* portable aluminium alloy box

Connection Type

WiFi (802.11n/20MHz/5G/450Mbps)
Probe built-in WiFi signal

Supporting Devices

Tablet, Smart phone, PC, etc

Working system

Android, IOS, Windows

Image save

jpg, avi, Mp4 and DICOM format

Buttons

3 programmable buttons

Lithium Battery

4200 mAh. 3.8V Support wireless charging

Working time: ≥3h, Standby time:≥12h

Battery life 5-8years

Dimensions

158(h)x70(w)x30(d) mm

Net weight 239g

Elements 192


Footprint R60

Boot-up

Platform dependent, generally <5 seconds



Offer great images to everyone





- * Training how to install the machine, training how to use the machine.
- * Spare parts replacement and other software/hardware trouble shooting
- * Offer OEM&ODM service for distributors
- * Offer SDK depends on different models
- * Free software and free upgrade the software
- * 18 Months Warranty, extendable up to 3 years

BEIJING KONTED MEDICAL TECHNOLOGY CO.,LTD

Tel:8610-60219113

Fax:8610-60219213

Service and Cooperation :info@konted.cn

Website:www.konted.cn Facebook: Beijing Konted

Youtube: Medical Konted

Toutube. Wedical Rollied

Address: Room 111,1F, Building 3, No. 27 Yongwang Road, Daxing Biological Pharmaceutical Industry Base, Daxing District, Beijing

102629, China.